

BeautiBond™

The Power Of One Unique Chemistry



BeautiBond™

BeautiBond is an innovative all-in-one 7th generation bonding system designed to increase efficiency and ensure a reliable bond between composite resins or GIOMER and natural tooth structure. The HEMA free formulation and incorporation of two unique adhesive monomers helps to save time, eliminate sensitivity and achieve predictable bonding with total peace of mind...

THE POWER OF ONE



UNIQUE DUAL ADHESIVE MONOMERS

Ensures a strong reliable consistent bond to both enamel and dentin

SINGLE APPLICATION

- Eliminates guess work
- Saves chair time

HEMA-FREE FORMULATION

- Eliminates post-operative sensitivity
- Provides stable & durable bond
- Minimized gingival blanching

PREDICTABILITY EVEN WITH A THIN BONDING LAYER

- Thin bubble-free bonding layer of 5µm
- Facilitates easy adaptation of restorative

CONVENIENCE

Available in Unit Dose or Bottle





SCIENCE BEHIND THE BOND

Two Powerful Adhesive Monomers

The dual adhesive monomers deliver almost equal bond strength to both enamel and dentin. The phosphonic acid monomer of **BeutiBond** is more stable than phosphoric acid monomer which enhances bonding to tooth enamel. The carboxylic acid monomer facilitates a strong bond to dentin.

SEM photographs courtesy of Dr. Toshimoto Yamada.



E – Enamel
D – Dentin
B – BeutiBond
CR – Composite Resin

Enamel Bonding Interface



Dentin Bonding Interface

Virtually No Post-Operative Sensitivity

The carboxylic acid monomer in **BeutiBond** helps demineralise the smear layer without exposing the dentinal tubules which virtually eliminates post-operative sensitivity.

Stable & Durable Bond

The HEMA-free formulation prevents aging caused by hydrolysis of the hybrid layer and adhesive layer to guarantee a stable bond with outstanding bonding durability.



Source – "Micro-tensile bond strength of one, two-step self-etching adhesives using novel thermal cycling method". Kawamoto Chiharu, Nagano Futami, Yasumoto Kei, Sano Hidehiko; Graduate School of Dental Medicine, Hokkaido University, Sapporo, Japan. "3rd International Congress on Adhesive Dentistry 2008".

Source – Internal data by SHOFU INC Japan.



Order Information



BEAUTIBOND SET [PN 1781]
BeautiBond 6ml bottle, Microtip Fine (with 1 handle) x 50 application tips, V-Dish x 25



BEAUTIBOND UNIT DOSE [PN 1782]
BeautiBond Unit Dose 0.1ml x 50, MicroBrush Fine x 50

Related Products



BEAUTIFIL II 6 COLOR SET
6 x Beautifil II Syringe, 4.5g – A2, A3, A3.5, B2, A3O, Inc, Dura-White Stone, Super-Snap Singles, OneGloss Trial Set & Paper Pad



BEAUTIFIL II COSMETIC SET
3 x Beautifil II Syringe, 4.5g – A2, A3, Inc, Beautifil Opaquer-LO, Beautifil Flow F02-A2, Dura-White Stone, Super-Snap Singles, OneGloss Trial Set & Paper Pad



BEAUTIFIL II REFILLS (4.5g)
Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, Inc, BW, A2O, A3O



BEAUTIFIL OPAQUER
(2g with needles tips x 5)
Shades: UO (Universal), LO (Light)



BEAUTIFIL FLOW – F02
F02 – Low Flow (2g with needles tips x 5)
Shades: A1, A2, A3, A3.5, A4, A3O, G




BEAUTIFIL FLOW – F10
F02 – High Flow (2g with needles tips x 5)
Shades: A1, A2, A3, A3.5, A4, A3O, A3T

CONVENIENCE & PREDICTABILITY

Thin Bonding Layer

The novel adhesive technology of **BeutiBond** creates a very thin bubble-free bonding layer of 5µm that helps to facilitate easy application of the first layer of the restorative material.



SEM photograph courtesy of Dr. Naotake Akimoto.

Simple Application


Etch, prime and bond in **just 30 seconds** with a single application and be assured of a reliable, strong bond to both enamel and dentin.



Apply **BeutiBond** (5 sec).
Leave undisturbed for 10sec.



Air-dry for 5sec to evaporate water and volatile components.
(Gently air-dry for 3sec then air-dry with high pressure for 2sec)



Light-Cure (Halogen: 10sec / LED: 5 sec).

Guaranteed Results


Photographs courtesy of Dr. Sushil Koirala.



Prepared cavity.



Application of **BeutiBond**.



Air-dried and light-cured.



Completed restoration with **BEAUTIFIL II**.





SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605 - 0983, Japan (www.shofu.co.jp)
SHOFU INC. Singapore Branch. 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
SHOFU DENTAL SUPPLIES (SHANGHAI) CO. LTD. No.645 Jiye Road, Sheshan Industry Park Songjiang Shanghai 201602, China (www.shofu.com.cn)
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)
SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8JU, UK (www.shofu.co.uk)

SHOFU is a registered trademark of SHOFU INC. All other trademarks and registered trademarks are the property of their respective holders. SHOFU INC. reserves the right to change specifications without notice.
© Copyright 2001 SHOFU INC. All rights reserved.